

Guangzhou Langwen Technology Co.,Ltd

 [+86-13829771029](tel:+86-13829771029)

 sales@loveingbeauty.com

 [RM508 ,Baoping Road, HaiZhu Area, GUANGZHOU CHINA](#)

Focus on beauty instrument OEM
and export source factory



PRODUCT BROCHURE

2024



Guangzhou Langwen Technology Co.,Ltd. more than 10 years original and powerful factory ,more than 10 years focusing on the development of beauty equipment, Pursuit of higher product quality, Making your business easier! is a professional beauty device company specializes in the development, production, sales and service of beauty equipment,

- * 755+808+1064nm Ice Cooling Diode Laser Series,
- * multi-function series 2 in 1 , 3 in 1, 6 in 1 series, Elight series, Q Switch ND YAG laser series,
- * Fat lossing :Cryolipolysis, RF series ,Cavitation and so on.

All the products are taking with CE/ISO certification ,suitable for world market , mainly used on North America South America , Europe ,more than 102 countrys .

The company provides long-term OEM processing, ODM customiztion services, and welcomes domestic and for-eign distribution agents and beauty agencies to customize. Welcome to visit the factory.

KH CATALOG PRODUCT CATALOG




360° CRYOLIPOLYSIS
FAT FREEZING MACHINE

01-06 Page



BUILD MUSCLE
REDUCED FAT

07-12 Page




GOLDEN EAGLE
ROLLER INSTRUMENT

13-18 Page



LOSE WEIGHT
SLIMMING DEVICE

19-23Page




BREAST ENHANCEMENT
BUTTOCK LIFTING
DEVICE

24-26 Page




TATOO REMOVAL
HAIR REMOVA
MACHINE

27-34 Page



FACIAL BEAUTY
TREATMENT

35-39 Page



RF MICRONEEDLE
BEAUTY MACHINE

40-41 Page



ARGUMENT

What is 360 Fat Freezing?

Fat freezing 360, also known as cryolipolysis, is a non-surgical way to reduce fat deposits in certain areas of the body. The procedure uses cold temperatures to kill fat cells and is designed for people who have stubborn fat that doesn't respond to diet and exercise.




Handle control method	Touch	Large cryo handle size:230*105mm	
Screen size	10.4	195*85mm 210*90mm 215*100mm	
Handle screen size	4.5	Small cryo handle size: 165*90mm	
Handle material	Silicom	140-70mm 155*80mm	
Handle temperature	-10°C/+45°C	Double chin handle size	80*40mm
Pressure output	0.08MPa	Cavitation head frequency	40KHz
Refrigeration method	Water cycle	Laser power	100mv/pcs
Handle face	3MHz	Laser wavelength	650nm
Handle frequency	5MHz	Number of laser lights	72PCS







NEW TECHNOLOGY --360 ANGLE CRYO



360° all-dimension surrounding cooling applicators
 The treating area is all-dimension cryolipolise Approach 100% frozen in full area
 A FREEZE CONTOUR PLUS treatment corresponds approximately 3 treatments with typical fat freezing devices Only 35 to 45 mins time each treatment needed instead of 60 mins . Save about half treatment time to treat more customers





- ➔ Fat Burn
- ➔ lifting
- ➔ tightening
- ➔ slimming
- ➔ shaping
- ➔ beauty



40K Cavitation

Body RF

Face RF

Lipo laser



WL-001-2

size:
60*45*65CM



WL-001-1

size:
60*45*65CM



WL-001-10

size:
76*43*63CM



WL-001-8

size:
58*48*70CM



WL-001-11

size:
56*50*55CM




WL-001-9

size:
66*53*67CM


HIEMT machine

By using HIEMT(high energy focused electromagnetic wave) technology to expand and contract muscles, limit training, and deep remodeling of muscle internal structure, That is, the growth of myofibrils(muscle enlargement) and the production of new protein chains and muscle fibers (muscle hyperplasia) to train and increase muscle density and volume.



How Does It Work?

E SCULPT is based on high-intensity focused electromagnetic energy. A single E SCULPT session causes thousands of powerful muscle contractions which are extremely important in improving the tone and strength of your muscles.



HIEMT 100% limit muscle contraction of the technology can lead to a large amount of fat decomposition, fatty acids from the decomposition of triglyceride outflow, a large amount of accumulation in adipocytes. Fat acidity is too high, fat cells apoptosis, in a few weeks by the body's normal metabolism out of the body.


RF: Heat uniformly for the treatment area, so that the subcutaneous fat quickly reaches the treatment temperature, achieving the dual effects of firming and dissolving fat.



WL-M10-2

Product Name: EMS Neo
Power: 2000W
Voltage: 110V-220V 50/60HZ
Treatment Handles:
Standard(4handles for body sculpting)
Optional(Sits stool for pelvic floor)
Frequency: 5hz~100hz
Dimension: 43X40X91CM (43kg)
Packing Size: 56X60X112 CM
weight: 64.30kg

SCHEMATIC DIAGRAM OF OBESE SKIN



HIGH ENERGY FOCUSED ELECTROMAGNETIC PRINCIPLE

Build muscle & Burn fat

Before treatment In treatment After treatment

Product name: EMS HIEMF Body Sculpt Machine
Feature: Weight Loss
Technology: HIEMF
Keywords: Muscle Sculpt Machine
Handles: 2 Handles
Treatment area: Abdomen Buttock
Voltage: 110V/220V 50-60Hz
Packing Size: 116X78X59CM
weight: 80.00kg



WL-M10-4




WL-M10-1

Product Name: EMS Neo
Power: 3000W
Voltage: 110V-230V 50/60HZ
Treatment Handles:
Standard(4handles for body sculpting)
Optional(Sits stool for pelvic floor)
Frequency: 5hz~100hz
Dimension: 54X65X104.5cm
Packing Size: 58X75X122.5cm



WL-M10-8

Power: 3000W-5000W
Magnetic vibration intensity: 7Tesla
Output power:300W-5000W
Output power:3-200HZ
Fuse:20A
Voltage: 110V/220V 50-60Hz
Dimension: 39X52X34cm
Packing Size: 46X64X79cm



BEFORE

Skin,fat and muscle form your overall body appearance.

DURING

The skin remains unaffected while the energy penetrates to fat and muscle layers.

AFTER

This results in substantial improvement of your appearance due to building muscle and burning fat.





WL-M10-3

Handles: 2 Handles or 4 Handles
Voltage: 110V/220V 50-60Hz
Handpiece: 4 pieces or 2 pieces
Pulse width: 300us
Function1: Weight Loss Beauty Machine
Function 2: Build Muscle+Fat Burning
Function 3: Body Fitness Machine
size: 44X55X65cm

FINE CARVING WITH ROLLER BELT G5

Inner roller Therapy is through the transmission of low frequency vibrations can generate a pulsed, rhythmic action on the tissues. The method is performed through the use of handpiece, selected according to the area of desired treatment. Time of application, frequency and pressure are three forces determining the treatment's intensity, which can be adopted to the clinical condition of a specific patient. The direction of rotation and the pressure used ensure that compression is transmitted to the tissues. Finally, it works to lift and firm, Cellulite Reduction, and lose weight.



SR-110



SR-111



- 1.6 technologies-80K Cavitation + RF + Vacuum + Micro-current + Heat + Magnetic vibration+lipolaser
- 2.80K cavitation-strong energy and better slimming effect.
- 3.3 vacuum RF handles-respectively for treatment on eye area, facial and body skin.
- 4.Large touch screen-easy and convenient to operate. 5.Multiple Functions-Face lifting, Body slimming and weight lose

Product name: Vacuum Cavitation Slimming Machine

Vacuum Cavitation Slimming Machine: Cavitation + RF + Vacuum +Micro-current + Heat + vibration+Lipolaser

Handpiece: 10

Cavitation Frequency: 80KHz

Vacuum RF handles: 3pcs (Small, Medium, Large)

Keywords: Cavitation Machine Ultrasonic

Function: Cellulite Reduce


Function : Body &face&eye Slimming Skin Lifting

Voltage: AC110V/220V 60-50Hz

Single package size: 117X50X52 cm

Single gross weight: 28.00 kg

BEFORE AND AFTER





GL-87T11

Product name:
Inner ball roller machine
Touch screen: 10.4 inch large LCD
Speed of large handle: 450rpm
Speed of small handle: 410rp
Output frequency::40-254H
Input voltage: AC110V/220V
Output power: 10-300W
Fuse: 5A
Package: Aluminium case
Package size: 56 X 67 X 76 cm
Gross weight: 40.00 kg



GL-87T10

Target Area: Hands, Body
Feature : Weight Loss, Skin Tightening, Cellulite Reduction
Handles: Facial*1 + Body*1
Speed of small handle: 410rpm
Speed of large handle: 450rpm
Output frequency: 40-254Hz
Output power: 100-1000W
Ambient temperature: 5° C~+40° C
size:55 X 40 X 65cm
weight:39.00kg

USE FOR THE BODY

LARGE HANDLE

Adjust speed button

LED show light pole

LED display screen

A roller device composed of 50 silicone balls in the larger handle



Honeycomb sphere arrangement

Clockwise rotation (Abbreviation: clockwise rotation)

Rotate counterclockwise (Abbreviation: reverse)



Target Area: Legs/Arms, Hands, Neck/Throat, Body, Face, Breast
Brand Name: Lislang
Application: For Commercial, Home and office use
Model Number: GL-87T12
TOUCH SCREEN: 5.6 inch LCD
SPEED OF LARGE HANDLE: 450rpm
SPEED OF SMALL HANDLE: 410rpm
OUTPUT FREQUENCY: 40-254Hz
INPUT VOLTAGE: AC 110V/220V
Single package size: 117X56X42 cm
Single gross weight: 42.00 kg



GL-87T12

MULTIFUNCTIONAL CAVITATION MACHINE

This Cavitation+RF equipment is the latest in modern technology. It combines cavitation with RF. Cavitation uses an effect to create air bubbles around the lipocyte membranes to break them apart. The lipocyte membranes will be absorbed and metabolized out of the body by the lymphatic system. Extremely effective to achieve weight loss. 6 MHz radio frequency is used to change the electric field inside the skin from positive to negative. The result will be a tighter, smooth and rejuvenated skin.



Model Number SR-114


Single package size: 38X32X35 cm

Single gross weight: 8.00 kg



40K fat explosion instrument:

crack fat, effectively consume heat, cell moisture, reduce fat cells, so as to achieve the effect of removing fat.



Model Number: SR-113
 Product Name: 5 in 1 cavitation rf
 Technology: Vacuum+Cavitation+ RF
 Handles: 4 Handles+8pcs Lazer pads
 Function: body slimming/shape/face lift
 Cavitation: 40K
 Color: white / pink
 Application: Beauty Salon
 Packing Size: 53X37X32cm
 Gross weight: 15.00KG





SR-100

Place of Origin: Guangdong, China
Model Number: SR-100
Type: Vacuum System
Keywords: kim 8 new rf vacuum
slimming machine
Style: Portable
Voltage: 110V/220V
Weight : 8.00kg
Packing size: 40X38X46cm
Function: 40k machine
Handle: 5 hands

Product name:Vortex Cavitation
RF Vacuum Machine
Technology: Vacuum+Cavitation+
RF+LED
Handles: 7 Easy Fast-connect
Style: Desktop
RF frequency: 1~5M Hz
Power: 200W
Voltage: 110V/220V 50-60Hz
Single package size: 58X40X35 cm
Single gross weight: 15.00 kg



SR-107

BREAST ENHANCEMENT

The specially designed vacuum pump can quickly improve size and shape of breast, without any aching and health threat, it does not involve silicon injection or pushing and filling a bra, it would make you have natural and pretty breast according to physiological period.



BE-MT05

Model Number: BE-MT05
Function 1: Breast Massage
Lifting Enhancement
Function2: Facial Lifting
Warranty: 1 year for main
machine,3 months for parts
Application: Beauty Salon+spa
Feature: Portable
Treatment area: Breast,Face
Package Size:37x36x38cm
Single gross weight: 6.00 kg





BE-MT01

Model Number: BE-MT01
Voltage: 110V/220V 50-60Hz
Technology: Vacuum+ Vibration+
Far Infrared
Negative Pressure: 60-70cmHg
Componet: 32 Cups + 3 Metal
Heads
Application: Beauty Salon + Spa
Single package size: 57X40X45cm
Single gross weight: 6.00 kg

Power: 15W
Vacuum: -73KPa
G.W.: 5KG
Package Size:
38x32x20cm(15x12.6x7.85inch)
Packing List
1 x main machine
3 pairs of large cups
3 pairs of cups for cupping
6 Pairs of gas cups for scrub



BE-MT06



HAIR REMOVAL LASER PICOSECOND MACHINES


808nm diode laser allows fast repetition rates up to 10Hz (10 pulses per second), fast hair removal for large area treatment.

755nm+808nm+1064nm 3 wavelength , more suitable for all skin color hair removal treatment.

1800w highest handle power, high laser energy density output up to 120J/cm²

Effectiveness. Maximum efficiency to achieve freezing and painless hair removal treatment.




1000w , 1200w , 1600w 1800w,2000W big power for free choice



We provide 5 treatment heads, one scanning head, one cutting head, three private heads
Scanning head is used for tender skin, shrinking pores, removing striae gravidarum and chloasma
Cutting head is used for removing sarcoma/nevus
Private head is used for vagina tighten, tender vulva, tender vulva



TX-M-12





TX-M-22

Model Number: TX-M-22
Spot size: HR: 16*57mm SR: 8*34mm
Wavelength: HR:690-1200nm
SR:560-1200nm
Energy: 2000W
Dimensions: 52.5X49X108cm
Net weight: 45.00kg

RF radio

Using unique RF technology and a contact cooling system, your skin can be significantly tightened and lifted, and subsequent gelatinization builds further tightens the dermal tissue and reduces wrinkles.



Lifting and tightening
 Promote collagen renewal



TX-M-12

Laser Type: Nd Yag Laser
Style: Stationary
Product name: picosecond machine
Wavelength: 1064nm/532nm
Pulse width: 500ps&800ps
Spot size: 2-10mm²
Pulse energy: 100-500mJ & 1064nm;
50-250mJ & 532nm
Pulse Rate: 1,2,5,10Hz



TX-M-11

power: 2000W
work: laser
wavelength: 1064nm; 755nm; 32nm
laser output frequency: 1-10hz
spot diameter: 2-10mm adjust
pulse: 800mj
size : 90X43X108cm
machine Net weight: 74.00kg

ND: YAG LASER


Using the Q-Switched ND:YAG model, the device uses a laser to instantly emit high energy, smashes pathological tissue pigments, and splits into particles. These particles are swallowed by macrophages and finally eliminated by the lymphatic system. The pigment gradually became lighter and eventually disappeared.



800w energy

Remove pigment spots

Equipped with three probes



Pulse hair removal

Slides with three wavelengths

particularly effective to hair follicle melanocytes without injury surrounding tissue. The laser light can be absorbed by hair shaft and hair follicles in the melanin, and converted into heat, thus increasing the hair follicle temperature. When the temperature raises high enough to irreversibly damage the hair follicle structure, which disappears after a period of nature physiological.

Gemstone crystal
 freezing point, painless





TX-M-21

Product name: Co2 Fractional Laser
Wavelength: 10600nm
Spot size: 0.02~0.05mm²
Working mode : Continuous operation
Distance: 0.1-2.6mm
Cooling system: air circulation
Energy: 1- 30W
Single package size: 82X107X57 cm
Single gross weight: 77.00 kg

Laser Type: Diode Laser
Laser Quantity: Germany Dilas bars
Laser wavelength: 8nm/755nm/1064nm
Spot Size: 15x15mm/15x30mm
Repetition Frequency: 1 ~ 20HZ
Pulse Width: 10-40ms
Energy Density: 1 ~ 30 J/cm
Laser Cooling Temperature: -5°C - 5 °C
Net weight 60.00kg
Package Size 68X52X70cm



TX-M-19



TX-M-10

Laser type: Picosecond laser
Wavelength: 532nm\1064nm\755nm\1320nm
Indicator light: 650nm infrared light
Maximum Energy: 1000MJ
Spot diameter: 2 -10mm
Working Frequency: 1 - 10HZ
Pulse width: 550-750ps
Power: 1000w
Single package size: 63X47X58 cm
Single gross weight: 25.00 kg

Laser Type: ND YAG Laser
Display: 10.4" TFT true color touch screen
Wavelength: 1064nm, 532nm, 755nm, 1320nm
Power supply voltage: AC220V; AC110V
Frequency: 1~10Hz
Pulse width: 6ns
Input power: 2000W
Packing size: 50X65X70cm
Weight: 27.00kg



TX-M-11



TX-M-17

Wavelength: 808nm/755nm/1064nm
Spot Size: 12*24mm²/12*36mm²
Pulse Duration: 10~400ms adjustable
Frequency: 1.0-10 Hz
Handle Power: 1200w/1600w/1800w
Machine Power: 2000w
Crystal temperature: -20~5° C
Single package size: 115X68X59 cm
Gross weight: 95.00kg

Laser type Semiconductor laser
wavelength 810nm
Spot area 10*12mm
Interface 8.4 inch touch screen
Input power 2200W (AC220V/ 10A, AC110V/20A)
Energy density 5 ~ 10J/cm²
Repetition rate 10Hz (fixed)



TX-M-20



TX-M-10

Power: 300W
Pulse Width: 20-950ms
Cooling System: Water+Wind+Tec
Frequency:1~10Hz
Voltage: AC 110~250V/50-60Hz
Spot Size: 15-30 mm
Single package size: 54X58X112 cm
Single gross weight: 50.00 kg

Out power: 500W
Wavelength: 808nm
Frequency:1-10HZ
Laser bar: 50W/Piece
Spot size: 10*10mm/15*15mm/15*20mm
Energy Density: 1~120J/cm²
Shots: 20 million shot
Laser wavelength: 808/755/1064
Single package size: 61X43X43 cm
Single gross weight: 55.00 kg



TX-M-14

FACIAL BEAUTY TREATMENT

Adsorbs the skin capsule, effectively cleans, dissolves cell wall waste, moisturizes and moisturizes.

Whitening skin and improving the skin. Stimulates skin tissue and hair follicles, resists skin aging and enhances skin elasticity.

Deep cleansing, while at the same time providing sufficient water, oxygen, and nutrient media to the skin.

Removes blackheads, shrinks pores, improves sagging skin and enhances skin clarity.



360 collagen gun

Promote collagen regeneration, improve fines around the eyes
promote body skin boby skin firmness and smoothness, reduce fat and shape



MZ-L2122

Collagen gun function

Improves favial contours and reduces wrinkles
Fat loss and body shaping
collagen regeneration


Facial head

Promote collagen regeneration, improve fine lines around the eyes, dark circles, bags under the eyes



Body head

It has a fat-melting mode, and the fat-melting can be turned on during the operation, and the effect of comprehensive firming and shaping can be achieved through the operation of the handle of the body.




BEFORE

Uneven skin appearance caused by loss of collagen due to natural aging process.



DURING

Heat energy disrupts old collagen fibers and initiates the process of new collagen formation.



AFTER

Treatment leads to visible skin improvement.



Model Number: AL-3004

Product Name: Multifunction Hydra
Dermabrasion Machine

Function: Moisturizing, Whitening,
Wrinkle Removal

Technology: RF+Photon+Vibration+Ion+Cool

Handle: 14 handles

Screen: 15 Inch Touch Screen

Single package size: 118X54X47 cm

Single gross weight: 60.00 kg

Functions:

Hydro Facial Probe Oxygen Sprayer
Skin Scrubber Hot Cold Hammer
Bio with Photon Ultrasonic Probe

Gross Weight: 3.60 KG

Packing Size: 38X26X16 CM



AL-3005



AL-3006

New Technology For Skin Rejuvenation

Package Listing:

- 1 x Main machine
- 1 x Hydro facial probe
- 1 x Ultrasonic probe
- 1 x Hot and cold hammer
- 1 x Skin scrubber
- 1 x Bio with photon probe
- 1 x Oxygen sprayer

Product Size : 38 x 31 x 35 cm


Gross Weight : 5.70kg

RF MICRONEEDLE BEAUTY MACHINE

What is gold radio frequency microcrystal

The gold radio frequency microcrystal is the ingenious combination of micro crystal and radio frequency. The two word of "gold" originates from the microcrystalline gold coating and the coating is also golden yellow. At the time of treatment, the doctor according to the crystal position on problems and treatment, adjust the depth of penetration and crocrystalline radiofrequency energy, then in the electronic control system, at the same time, dozens of insulating ceramics penetrate the skin quickly, from the micro crystal tip radiofrequency energy, then quickly exit, so the cycle until the treatment is finished.

9 steps to get watery clear skin





SR-RU126

Screen: 8.0 inch LCD touch screen
Feature: Face Lift, Skin Rejuvenation, Wrinkle Remover
Needles: 4pcs (10pin, 25pin, 64pin, Nanocrystals)
Voltage: AC110V/220V 50-60Hz
Output Power: 10-200W
RF frequency: 2-4MHZ
needle depth: 0.2-3.5mm(0.1 step)
Package size: 49X48X80cm
G.W: 15KGS

Function: Scar Removal
Keywords: Microneedling
Screen: 8.4 Inch Color Touch Screen
Voltage: 110V/220V 50-60Hz
Power: 200W
Technology: microneedle RF+ Fractional RF
RF Frequency: 5Mhz
Single package size: 118X50X60 cm
Single gross weight: 65.00 kg



Muscle Happiness Chair

Non-invasive private beauty equipment

<u>energy: 90J/cm²</u>	<u>depth: lamina propria 5mm</u>
<u>time: single shot for 8 seconds</u>	<u>Intrinsic temperature: 45°C-50°C</u>
<u>mucosal temperature: 38°C-41°C</u>	<u>size: 83X80X91.cm</u>
<u>weight: 109.00kg</u>	<u>logo: customized Logo available</u>

